

Pressure Regulator Habicht with Tube Flowmeter



Manufacturer: Greggersen

Product 910003

Number:

Description

O2 Pressure Regulator Unit

HABICHT

Suction and Insufflation Unit,

Cylinder Connection

Combi-unit with cylinder connection, pressure regulator with high-pressure cylinder connection according to DIN 477, preset to 500 kPa (5 bar), with contents gauge, safety valve, ejector with vacuum gauge, flowmeter with fine dosing valve for regulating the flow rate;

Scope of delivery without overflow jar or humidifier bottle

Gas type: Oxygen / O2

Connection: High-pressure cylinder connection according to DIN 477, No. 9

Outlet: Tubing Nipple / Vacuum Nipple

Output: 0...15 l/min

adjustable between 0 and -85 kPa (0 to -0.85 bar)

Consisting of:

1 x Combi-Unit