


Kolibri, Rail O2 Standard NIST ISO32



Manufacturer: Greggersen

Product 902164

Number:

Description

Kolibri, O₂, 0...15 l/min

RAIL-MOUNTED UNIT NIST

Medical Compact Flowmeter according to DIN EN ISO 15002;

Flow rate is adjusted via handwheel in 10 flow settings;

Gas inlet NIST integrated into rail clamp;

Outlet 9/16-18 UNF and tubing nipple

Handwheel in ISO32 - white for O₂

Gas type: Oxygen / O₂

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF - 18 UNF / tubing nipple

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
