

Pirol 2.0 -90,VAC,PROBE, MEDAP



Manufacturer: Greggersen

Product 904613

Number:

Description

Vacuum Regulator Pirol -90, Plug-in Unit MEDAP

Medical vacuum regulator for vacuum dosing.

With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel.

The adjustment mechanism enables infinitely variable, precise and intuitive regulation of the vacuum. The rotatable vacuum gauge allows the user to read the set value from an optimal viewing angle.

The regulator unit meets the requirements of DIN EN ISO 10079-3 (high vacuum / high flow).

Material:

Body: aluminium, anodised

Handwheel: plastic

Probe: brass, chrome-plated

Technical data:

Gas type: vacuum / VAC

Inlet: MEDAP VAC Probe

Control range: 0 to -90 kPa $\pm$ 10%

Flow rate: >20 l freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum
-----------------	------------
