


## Probe with Hose Barb N2, AGA



**Manufacturer:** Greggersen

**Product** 505166

**Number:**

### Description

Probe with Hose Barb for Nitrogen  
with 7 mm Hose Barb and Non-Return Valve

Gas Type: N2

Coding: N2

Inlet: Probe to SS 87 524 30 (AGA)

Outlet: 7 mm Hose Barb

### Properties

|                 |             |
|-----------------|-------------|
| <b>Gas Type</b> | N2-Nitrogen |
|-----------------|-------------|