

Double Flowmeter, MEDAP, O₂, 0-15l/min



Manufacturer: Greggersen

Product 901577

Number:

Description

Double Flowmeter, O₂, 0...15 l/min, MEDAP

-Plug-in Unit-

Medical double flowmeter, pressure-compensated according to DIN EN ISO 15002, with fine control valve for dosing the flow, for coupling into a MEDAP terminal unit

Gas type: Oxygen / O₂

Connection: MEDAP Probe

Outlet: 9/16-18 UNF

Flow rate: each 0...15 l/min

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type

O₂-Oxygen