

Hose, N2O, 5m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900471

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled probe with reverse flow limiter,
hose and NIST device connection

Specification:

Gas type: Nitrous oxide / N2O

Inlet: Probe according to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type

N2O-Nitrous Oxide