

Nebuliser Plug-in Valve with Clamp AIR/O₂



Manufacturer: Greggersen

Product 900865

Number:

Description

Fitis Nebuliser Plug-in Valve with Clamp AIR/O₂

Plug-in Unit

Nebuliser plug-in valve, fixed setting with side clamp and

Combi-Probe, for connection to a terminal unit

Gas type: Compressed Air / Oxygen

Connection: Probe to DIN 13260 Part 2

Outlet: 9/16-18 UNF/ Tubing Nipple

Flow Rate: 5 l/min ± 20%