

## Falke with Flow Standard O2 PIN 9/16



**Manufacturer:** Greggersen

**Product** 500027

**Number:**

### Description

Falke with Flow, Standard,  
Compact pressure regulator with high-pressure connection  
acc. to CGA V-1, contents gauge, step-adjustable  
flow rate.

Gas type: Oxygen / O<sub>2</sub>

Connection: HP cylinder connection acc. to CGA V-1; Connection: 870 (pin-index)

Outlet: 9/16-18 UNF (tubing nipple), bottom

Capacity: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

### Properties

**Gas Type** O<sub>2</sub>-Oxygen