

Hose, O2, 1.5m, NF - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 902526

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and NIST device connection

Specification:

Gas type: Oxygen / O₂

Inlet: Probe according to NF S 90-116

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type O₂-Oxygen