


Forano, Front Part, BS, O2, ISO



Manufacturer: Greggersen

Product 905160

Number:

Description

FORANO Terminal Unit - Front Part BS

Gas type: Oxygen / O2

Labelling: oxygen / O2

Gas-specific socket assembly for mounting
on FORANO base part.

Socket assembly according to BS 5682

Release via release ring.

Colour coding:

Socket assembly coloured according to DIN EN ISO 9170-1 (ISO32).

Release ring labelling: printed

Cartridge incl. sealing elements

Technical data:

Gas type: Oxygen / O₂

Labelling: oxygen / O₂

Colour coding: according to DIN EN ISO 9170-1 (ISO32)

Front diameter: Ø 42mm

Scope of delivery:

- Gas-specific socket assembly
- Release ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type

O₂-Oxygen
