

Automatic Manifold MediControl MC2050T



Manufacturer: Greggersen

Product 326051

Number:

Description

Automatic Manifold MediControl 2050T

Flow Rate: 50 m³/h

Inlet Pressure: max. 200 bar

Outlet Pressure: 5 bar +/- 0.5 bar

Electr. Connection: 230 V AC

System Language: German

Mode: Tank and 2 cylinder banks

Gas Type: ++++++

Consisting of:

1 x MC2050

1 x Pipe Set MC2050

1 x Pressure Regulator Tank Panel (327800)

Microprocessor-controlled gas supply system for 2-sided cylinder supply, design according to DIN EN ISO

7396-1, LCD display indicating operating status and service messages in plain text, forwarding of messages via potential-free contacts, sensory pressure monitoring, 2-stage pressure reduction, 2nd stage configured redundantly, pneum. priority switching in the event of power failure, with pressure gauges and emergency supply facility (NIST), completely mounted on mounting plate incl. cover; Additionally: Pressure regulator panel for tank supply with shut-off valve, pressure sensor, emergency supply facility (NIST)