

Compressed-Air Vacuum Regulator Specht -90 Probe

ANTI



Manufacturer: Greggersen

Product 500244

Number:

Description

Vacuum Regulator Specht-90, DIN Plug-in Unit, non-magnetic

Medical compressed air suction unit for vacuum generation, with Vacuum Gauge, fine dosing and quick-action shut-off valve, silencer and filter holder.

Technical data:

Gas type: Compressed air / AIR

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------