


Falke without Flow O2 DIN +1DIN



Manufacturer: Greggersen

Product 910382

Number:

Description

Falke, without flow with 1 quick-action coupling,
Compact pressure regulator with contents gauge;
with one quick-action coupling at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Quick-action coupling according to DIN 13260-2, bottom

Output: 450 kPa (4.5 bar)

Properties

Gas Type O2-Oxygen