

Kolibri 2 Probe Carba O2 Standard+



Manufacturer: Greggersen

Product 906457

Number:

Description

KOLIBRI 2, O₂, 0...30 l/min

- Standard+ CARBA PLUG-IN UNIT -

Flow settings:

0 / 2 / 3 / 4 / 6 / 8 / 12 / 15 / 20 / 25 / 30 l/min

Compact Flowmeter KOLIBRI 2

with permanently attached orifices with selection device

Medical compact flowmeter for connection

to a terminal unit, with handwheel - incrementally

adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and legibility

of the set flow setting from various viewing directions

and angles possible.

Frontal operating concept -> increased safety when
setting the correct flow rate

Gas type: Oxygen / O₂

Connection: CARBA Probe

Outlet: 9/16-18 UNF

Flow rate: 0...30 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	CARBA
Variant	Standard+