

Kolibri, DIN Probe, AIR Standard+ long



Manufacturer: Greggersen

Product 901486

Number:

Description

Kolibri, AIR, 0...30 l/min, Standard+, DIN

- PLUG-IN UNIT, long -

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
gradually adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas type: Compressed air / AIR

Connection: Probe according to DIN 13260 Part 2

Total length: Probe front part to centre of flow outlet 108 mm

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...30 l/min

Flow settings: 0/ 0.5/ 1/ 2/ 4/ 6/ 8/ 10/ 12/ 15/ 30+ l/min

Properties

Gas Type	AIR Compressed Air
----------	--------------------
