


Forano, Front Part, AGA, VAC, ISO



Manufacturer: Greggersen

Product 905142

Number:

Description

Terminal Unit FORANO - Front Part AGA

Gas type: Vacuum / VAC

Labelling: vacuum / VAC

Gas-specific socket assembly for mounting
on base part FORANO. Design complies with
DIN EN ISO 9170-1;

Socket assembly according to SS 87 524 30

Release via release ring.

Colour coding according to ISO 32.

Release ring labelling: printed

Cartridge incl. sealing elements

Technical data:

Gas type: Vacuum / VAC

Labelling: vacuum / VAC

Colour coding: ISO 32

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed labelling
- cartridge with sealing elements

Properties

Gas Type

VAC-Vacuum
