

Flowmeter, AGA, O₂, 0-15l/min



Manufacturer: Greggersen

Product 901655

Number:

Description

Flowmeter, standard

Medical flowmeter, pressure-compensated according to DIN EN ISO 15002, with fine control valve for metering the flow, for coupling into a terminal unit

PLUG-IN UNIT

Gas type: Oxygen

Connection: Probe according to SS 87 524 30 (AGA)

Outlet: 9/16-18 UNF

Flow Rate: 0...15 l/min

Properties

Gas Type

O₂-Oxygen