

Forano, Front Part, DIN, AIR, Colour-Neutral



Manufacturer: Greggersen

Product 905111

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral - Version III from 08/2019

Gas type: Compressed Air / AIR

Labelling: medical air / AIR

Gas-specific Socket Assembly for mounting
on Base Part FORANO.

Socket Assembly in accordance with DIN 13260 Part 2
with parking and operating position.

Release via Release Ring.

Colour Coding:

Colour-Neutral (black)

Release Ring Labelling: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas type: Compressed Air / AIR

Labelling: medical air / AIR

Colour Coding: Colour-Neutral (black)

Front Diameter: Ø 42mm

Scope of Delivery:

- gas-specific Socket Assembly
- Release Ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------
