


Kolibri, Carba Probe, O2 Standard



Manufacturer: Greggersen

Product 901585

Number:

Description

Kolibri, O2, Standard, 0...15 l/min, CARBA

PLUG-IN UNIT

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
stepwise adjustable flow rate, outlet 9/16-18 UNF and tubing nipple

Gas Type: Oxygen

Connection: Carba Med Probe

Outlet: 9/16-18 UNF / Tubing Nipple

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O2-Oxygen