

Falke w. Flow Stand. O2 NF Tubing Nipple BE-BU



Manufacturer: Greggersen

Product 910332

Number:

Description

Falke with Flow, Standard,
Compact Pressure Regulator with HP Cylinder Connection
according to NF E29-650, contents gauge, stepwise adjustable
flow rate.

Gas type: Oxygen / O₂

Connection: HP Cylinder Connection according to NF E29-650; Connection: F-627 (Bene)

Outlet: Tubing Nipple bottom

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O₂-Oxygen