


AGSS Forano AIR Rail Type1L



Manufacturer: Greggersen

Product 902075

Number:

Description

AGSS FORANO Typ1L
for rail mounting

according to DIN EN ISO 9170-2, with mounting bracket,
visual indicator for operational monitoring, all-metal
design, for operation with compressed air 5 bar,
external ejector with mounting bracket and
outlet pipe 15x1, incl. 500 mm hose each 4/6/15 mm

Connection: Cu pipe 8x1

Outlet: Cu pipe 15x1 (exhaust air)

6 mm coupling (drive gas)

Properties

Gas Type AIR Compressed Air