

Hose, AIR, 1.5m, DIN - NIST, ISO



Manufacturer: Greggersen

Product 902414

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow restrictor,
hose and NIST device connection

Specification:

Gas type: Compressed Air / AIR

Inlet: Probe according to DIN 13260 Part 2

Hose length: 1.5 m

Colour coding: ISO 32 (black/white)

Outlet: gas-specific NIST screw connection

Properties

Gas Type AIR Compressed Air