

Falke with Flow Neo O2 AGA Tubing Nipple



Manufacturer: Greggersen

Product 910336

Number:

Description

Falke with Flow, NEO,
Compact pressure regulator with high-pressure cylinder connection
acc. to SS 367615, contents gauge, incrementally adjustable
flow rate.

Gas type: Oxygen / O₂

Connection: High-pressure cylinder connection acc. to SS 367615; No. E1 (W21.8)

Outlet: Tubing nipple at bottom

Flow rate: 0...1 l/min.

Flow settings: 0/ 0.1/ 0.15/ 0.2/ 0.25/ 0.3/ 0.4/ 0.5/ 0.6/ 0.8/ 1.0 l/min

Properties

Gas Type O₂-Oxygen