


Falke w. Flow Stand+ O2 DIN 9/16



Manufacturer: Greggersen

Product 910351

Number:

Description

Falke with Flow, Standard+

Compact pressure regulator with contents gauge;
incrementally adjustable flow rate.

Gas type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (Tubing Nipple), bottom

Flow Rate: 0...30 l/min

Flow Settings: 0/ 0.5/ 1/ 2/ 4/ 6/ 8/ 10/ 12/ 15/ 30 l/min

Properties

Gas Type O₂-Oxygen