

Falke with Flow Paediatric CO2 DIN 9/16



Manufacturer: Greggersen

Product 500152

Number:

Description

Falke with Flow, Paediatric,
Compact pressure regulator with contents gauge;
incrementally adjustable flow rate.

Gas type: Carbon dioxide / CO2

Connection: Cylinder connection for gas cylinders according to DIN 477, no. 6

Outlet: 9/16-18 UNF (Tubing Nipple), bottom

Flow rate: 0...5 l/min

Flow settings: 0/ 0.1/ 0.3/ 0.5/ 0.6/ 0.8/ 1/ 2/ 3/ 4/ 5 l/min

Properties

Gas Type CO2-Carbon Dioxide