

Angled Probe AIR Hose Barb 7mm, DIN,ISO V4A



Manufacturer: Greggersen

Product 505151

Number:

Description

Angled Probe for Compressed Air
with 7 mm hose barb and non-return valve
Probe front part made of 1.4404 (stainless steel)

Gas type: Compressed Air / AIR
Connection: Probe according to DIN 13260 Part 2
Outlet: 7 mm Hose Barb
Colour coding: according to ISO 9170-1 (black/white)

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------