

Falke w. Flow Stand. AIR DIN Tubing Nipple +1DIN



Manufacturer: Greggersen

Product 500198

Number:

Description

Falke with Flow, Standard
with 1 quick coupling,

Compact pressure regulator with contents gauge;
incrementally adjustable flow rate,
DIN quick coupling at a pressure of 450 kPa.

Gas Type: Compressed Air / AIR

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 13

Outlet: Tubing Nipple / 1 quick coupling according to DIN 13260 Part 2

Flow Rate: 0...15 l/min Tubing Nipple - 120 l/min quick coupling

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type AIR Compressed Air