


Ventus 22 mm Valve Insert VAC



Manufacturer: Greggersen

Product 900782

Number:

Description

Valve insert for Ventus 22 mm

Gas type : VAC

Design complies with DIN EN ISO 7396-1

Valve assembly consisting of:

- Shut-off valve DN 20, G3/4"
- Copper pipe 22 mm inlet / outlet
- Pressure gauge -1...0 bar
- Pressure sensor -1...6
- NIST emergency feed point
- Physical separation point
- Gas type labelling

Properties

Fulfilled Standard	DIN EN ISO 7396-1
Gas Type	VAC-Vacuum
NIST Emergency Inlet Point	1 NIST Emergency Inlet Point