

Kolibri, Probe, DIN, AIR OFF/6



Manufacturer: Greggersen

Product 901402

Number:

Description

Kolibri, for Compressed Air

-Plug-in Unit-

Medical Compact Flowmeter according to DIN EN ISO 15002,
for connection to a Terminal Unit, with handwheel
step-adjustable Flow Rate, outlet 9/16-18 UNF
and Tubing Nipple

Gas type: Compressed Air

Connection: Probe according to DIN 13260 Part 2

Overall length: 88 mm (Probe Front Part to centre of Tubing Nipple outlet)

Outlet: 9/16-18 UNF / Tubing Nipple

Flow Rate: 0 / 6 l/min

Flow Settings: 0 + 6

Properties

Gas Type	AIR Compressed Air
Variant	Air Off/6