

## Falke w. Flow Stand. CO2 DIN 9/16



**Manufacturer:** Greggersen

**Product** 910390

**Number:**

### Description

Falke with Flow, Standard,  
Compact Pressure Regulator with contents gauge;  
incrementally adjustable flow rate.

Gas type: Carbon dioxide / CO2

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 6

Outlet: 9/16-18 UNF (Tubing Nipple) bottom

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

### Properties

**Gas Type** CO2-Carbon Dioxide