


Falke without Flow O2 DIN +1NF



Manufacturer:	Greggersen
Product	500353
Number:	

Description

Falke, O2, without flow
 with 1 NF quick-release coupling
 Compact pressure regulator with contents gauge
 and one quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: Quick-release coupling according to NF S 90-116, bottom

Output: 450 kPa

Properties

Gas Type	O2-Oxygen
-----------------	-----------