

Forano Base Part straight pipe, DIN, AGA, AIR8



Manufacturer: Greggersen

Product 903204

Number:

Description

Terminal Unit Forano Base Part DIN/AGA/NF

Gas-specific base part for installation in
supply rails or ducts.

Base part with spacer washer and mounting bracket;
brazed copper pipe, integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with Cu pipe
- spacer washer
- mounting bracket

Technical data:

Gas type: Compressed air 8 bar / AIR-800

Installation type: Media supply rail

Inlet: Cu pipe 8x1 mm straight, top

Properties

Gas Type

AIR-800-Compressed Air 8bar