

Hose, VAC, 3m, DIN - NIST, ISO



Manufacturer: Greggersen

Product 902417

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe,

Hose and device connection NIST

Specification:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Hose length: 3 m

Colour coding: ISO 32 (yellow)

Outlet: gas-specific screw connection NIST

Properties

Gas Type

VAC-Vacuum