


Kolibri, Carba Probe, O2 Standard+



Manufacturer: Greggersen

Product 901586

Number:

Description

Kolibri, Standard+

Medical compact flowmeter according to DIN EN ISO 15002,
for connection to a terminal unit, with handwheel
incrementally adjustable flow rate, outlet 9/16-18 UNF and
tubing nipple

PLUG-IN UNIT

Gas type: Oxygen

Connection: Carba Med probe

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...30+ l/min

Properties

Gas Type

O2-Oxygen
