

Forano, Front Part, DIN, CO2, Colour-Neutral



Manufacturer: Greggersen

Product 905115

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral - Version III from 08/2019

Gas Type: Carbon Dioxide / CO2

Labelling: carbon dioxide / CO2

Gas-specific Socket Assembly for installation
on Base Part FORANO.

Socket Assembly according to DIN 13260 Part 2
with parking and operating position.

Release via Release Ring.

Colour coding:

Colour-Neutral (black)

Labelling Release Ring: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas Type: Carbon Dioxide / CO₂

Labelling: carbon dioxide / CO₂

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of Delivery:

- gas-specific Socket Assembly
- Release Ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type	CO ₂ -Carbon Dioxide
-----------------	---------------------------------
