


Forano, Front Part, DIN,O2 93, Colour-Neutral



Manufacturer: Greggersen

Product 905118

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral

Gas Type: Oxygen 93 / O2 93

Labelling: oxygen 93 / O2 93

Gas-specific socket assembly for mounting
on Base Part FORANO.

Socket assembly according to DIN 13260 Part 2
with parking and operating position.

Release via release ring.

Colour coding:

Colour-Neutral (black)

Release ring inscription: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas Type: Oxygen 93 / O2 93

Labelling: oxygen 93 / O2 93

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of Delivery:

- Gas-specific socket assembly
- Release ring incl. printed inscription
- Cartridge with sealing elements

Properties

Gas Type	O293-Oxygen 93
-----------------	----------------
