

Vacuum Regulator Specht-90, AIR without Probe



Manufacturer: Greggersen

Product 500258

Number:

Description

Compressed-Air Vacuum Regulator Specht-90, connection G 1/4" male

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine adjustment and quick-action shut-off valve, silencer and filter holder.

Technical data:

Gas type: Compressed air / AIR

Inlet: G 1/4" male

Control range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type AIR Compressed Air