

Varioport Spatz -90



Manufacturer:	Greggersen
----------------------	------------

Product	902120
----------------	--------

Number:	
----------------	--

Description

Portable Suction Unit Varioport with Vacuum Regulator Spatz -90

Medical, portable suction unit for vacuum
for accommodating 2 secretion collection canisters.

The canister holder is universal. Common canisters in
sizes 0.5 to 3.0 L from various manufacturers can be used
(e.g. Serres, Abbott, Medela, Ardo, Cardinal, and many more).

The carrying frame features a sturdy handle and a bracket for
a 25x10 mm equipment rail. The vacuum regulator unit is located
at the front, allowing the operator quick and direct access.

The gas inlet as a NIST coupling is located at the rear of the
unit, ensuring that the supply hose remains outside the working area.

Scope of Delivery:

- Vacuum Regulator with fine-adjustment and quick-action shut-off valve

- Vacuum Gauge
- Carrying frame, coated
- 2 highly flexible retaining cuffs

Technical Data:

Gas Type: Vacuum / VAC

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF with vacuum tubing nipple

Control Range: -20 to -85 kPa

Weight: approx. 3 kg (without canisters)