

Kolibri, DIN Probe, O2 NEO



Manufacturer: Greggersen

Product 901405

Number:

Description

Kolibri, O2, 0...1 l/min

- NEO / PLUG-IN UNIT -

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
stepwise adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas type: Oxygen / O2

Connection: Probe according to DIN 13260 Part 2

Overall length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18UNF / tubing nipple

Flow rate: 0...1 l/min

Flow settings: 0 / 0.1 / 0.15 / 0.2 / 0.25 / 0.3 / 0.4 / 0.5 / 0.6 / 0.8 / 1.0 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
