

Kolibri 2 Probe DIN O2 Standard long



Manufacturer: Greggersen

Product 906404

Number:

Description

KOLIBRI 2, O₂, 0...15 l/min

- Standard DIN PLUG-IN UNIT -

Long version

Flow settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with permanently attached orifices with selection device

Total length "long"

Probe front part to centre of flow outlet

(distance from terminal unit to tubing nipple):

75 mm (intermediate piece 60 mm)

Medical compact flowmeter for coupling

into a terminal unit, with handwheel - incrementally adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability of the set flow setting possible from various viewing directions and angles.

Frontal operating concept -> greater safety when setting the correct flow rate

Gas type: Oxygen / O2
Connection: Probe according to DIN 13260 Part 2
Outlet: 9/16-18 UNF
Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:
900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O2-Oxygen
Inlet Standard	DIN
Variant	Standard Long