

Nebuliser Plug-in Valve with Socket O2



Manufacturer: Greggersen

Product 900775

Number:

Description

Fitis Nebuliser Plug-in Valve with Mount O2

Plug-in Unit

Nebuliser plug-in valve, pre-set with mount (top) and
Combi-Probe, for connection to a terminal unit.

Technical Data:

Gas type: Oxygen

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow rate: 5 l/min $\pm$ 20%

Properties

Gas Type

O2-Oxygen