

Falke w. Flow Stand+ O2 DIN Tubing Nipple +2DIN



Manufacturer: Greggersen

Product 910209

Number:

Description

Falke with Flow, Standard+

with 2 quick couplings,

Compact Pressure Regulator with contents pressure gauge;

step-adjustable Flow Rate.

DIN quick coupling at a pressure of 450 kPa.

Gas Type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing Nipple / 2 quick couplings according to DIN 13260 Part 2

Flow Rate: 0...30 l/min Tubing Nipple - quick coupling max. 120 l/min

Flow Settings: 0 / 0.5 / 1 / 2 / 4 / 6 / 8 / 10 / 12 / 15 / 30+ l/min

Properties

Gas Type

O₂-Oxygen