


Ventus 22mm 5-way Front Part+Base Part



Manufacturer: Greggersen

Product 900881

Number:

Description

Area Control Unit VENTUS 22mm

5-way, Front Part + Base Part

Gas type : *** Please specify ***

Language Instructions for Use: German, English, Dutch, Danish, Hungarian (please delete as applicable!)

Design complies with DIN EN ISO 7396-1 and HTM 02-01.

Front part assembly consisting of:

- Sheet steel door, powder-coated
- Side-hinged opening
- Alarm unit 2-way
- Magnetic foil to cover the warning
- Plaster compensation up to 25mm
- Power supply unit 230 V AC

Base part assembly consisting of:

- Built-in housing for flush mounting

The alarm unit monitors:

- Mains pressure too high
- Mains pressure too low

Dimensions

Inlet & Outlet: 22 mm CU pipe

Front part (WxHxD): 630 x 530 x 9 mm

Base part (WxHxD): 570 x 470 x 77 mm

Alarm limits of the warning:

O₂, AIR, N₂O, CO₂: 500 kPa ± 20%

VAC: -40 kPa

AIR-800: 800 kPa ± 20%