

Double Flowmeter, NIST, O₂, 0-15l/min



Manufacturer: Greggersen

Product 902620

Number:

Description

Double Flowmeter, O₂, 0...15 l/min

RAIL-MOUNTED UNIT NIST

Medical double flowmeter, pressure-compensated according to DIN EN ISO 15002, with fine control valve for regulating the flow; gas inlet NIST integrated into rail clamp

Gas type: Oxygen / O₂

Inlet: NIST according to DIN EN ISO 18082

Outlet: 2 x 9/16- 18 UNF

Flow rate: each 0...15 l/min

Properties

Gas Type

O2-Oxygen
