

Falke with Flow Stand. O2 BS Tubing Nipple Bull



Manufacturer: Greggersen

Product 910324

Number:

Description

Falke with Flow, Standard,
Compact pressure regulator with high-pressure cylinder connection
to BS 341-3, contents gauge, incrementally adjustable
flow rate.

Gas type: Oxygen / O₂

Connection: HP cylinder connection to BS 341-3; Connection: 3 (Bene)

Outlet: Tubing nipple bottom

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O₂-Oxygen