


Ventus 22 mm Valve Insert O2



Manufacturer: Greggersen

Product 900780

Number:

Description

Valve insert for Ventus 22 mm

Gas type: O2

Design complies with DIN EN ISO 7396-1

Valve assembly consisting of:

- Shut-off valve DN 20, G3/4"
- Cu pipe 22 mm inlet / outlet
- Pressure gauge 0...16 bar
- Pressure sensor -1...6
- NIST emergency inlet point
- Physical disconnection point
- Gas type labelling

Properties

Fulfilled Standard	DIN EN ISO 7396-1
---------------------------	-------------------

Gas Type	O2-Oxygen
-----------------	-----------

NIST Emergency Inlet Point	1 NIST Emergency Inlet Point
-----------------------------------	------------------------------