

Forano, Front Part, DIN, VAC, Colour-Neutral



Manufacturer: Greggersen

Product 905112

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral - Version III from 08/2019

Gas Type: Vacuum / VAC

Inscription: vacuum / VAC

Gas-specific Socket Assembly for mounting
on Base Part FORANO.

Socket Assembly according to DIN 13260 Part 2
with park and operating position.

Release via Release Ring.

Colour Coding:

Colour-Neutral (black)

Release Ring Inscription: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas Type: Vacuum / VAC

Inscription: vacuum / VAC

Colour Coding: Colour-Neutral (black)

Front Diameter: Ø 42mm

Scope of Delivery:

- gas-specific Socket Assembly
- Release Ring incl. printed inscription
- Cartridge with sealing elements

Properties

Gas Type	VAC-Vacuum
-----------------	------------
