

Kolibri, Probe, DIN, O2 Paediatric



Manufacturer: Greggersen

Product 901404

Number:

Description

Kolibri, O2, 0...5 l/min

- Paediatric / PLUG-IN UNIT -

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
incrementally adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Gas type: Oxygen / O2

Connection: Probe according to DIN 13260 Part 2

Total length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...5 l/min

Flow settings: 0 / 0.1 / 0.3 / 0.5 / 0.6 / 0.8 / 1 / 2 / 3 / 4 / 5 l/min

Properties

Gas Type	O2-Oxygen
Variant	Paediatric