

Falke with Flow Stand. O2 DIN 9/16 +2DIN



Manufacturer: Greggersen

Product 910362

Number:

Description

Falke with flow, standard with 2 quick-release couplings,
Compact pressure regulator with contents gauge;
incrementally adjustable flow rate,
DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (Tubing Nipple) / 2 quick-release couplings according to DIN 13260-2

Flow rate: 0...15 l/min 9/16-18 UNF; quick-release coupling max. 120 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type

O2-Oxygen