

Angled Probe N2O Hose Barb 7mm, DIN ISO V4A



Manufacturer: Greggersen

Product 505153

Number:

Description

Angled probe for nitrous oxide
with 7 mm hose barb and non-return valve
Probe front part made of 1.4404 (stainless steel)
Gas type: Nitrous oxide / N2O
Connection: Probe according to DIN 13260 Part 2
Outlet: 7 mm Hose Barb
Colour coding: according to ISO 9170-1 (blue)

Properties

Gas Type	N2O-Nitrous Oxide
-----------------	-------------------